

Insulin

For the Quantitative Determination of Human Insulin in Serum or Plasma

Cat. No. KAI-071

INTENDED USE

For the quantitative determination of human insulin in serum and plasma by immunoturbidimetric assay. FOR *IN VITRO* DIAGNOSTIC USE.

SUMMARY

Insulin is a peptide hormone with approximate molecular weight of 5,800 daltons. Secreted from β cells of the islet of Langerhans in the pancreas, insulin acts to reduce the blood sugar level. Since the blood insulin level reflects the function of β cells, insulin has been widely used as an important diagnostic tool for diabetes mellitus. The **K-ASSAY®** Insulin test is a highly specific assay for insulin in serum or plasma.

PRINCIPLE OF TEST

Latex particles coated with antibody specific to human insulin form immune complexes in the presence of insulin from the sample. The immune complexes cause an increase in light scattering that is proportional to the concentration of insulin in the serum or plasma sample. The light scattering is measured by reading turbidity at 600 nm primary, 800 nm secondary. The sample insulin concentration is determined versus insulin calibrators of known concentrations.

KIT COMPOSITION

Reagents (Liquid Stable)

R1: Buffer Reagent, pH 8.2
Tris(hydroxymethyl)aminomethane (100 mM)

R2: Latex Suspension
Anti-human insulin mouse monoclonal antibody (~1 mg/mL)

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use Only. R only.

Not to be used internally in humans or animals. Normal precautions exercised in handling laboratory reagents should be followed.

Do not mix or use reagents from one test kit with those from a different lot number.

Do not use reagents past their expiration date stated on each reagent container label.

Do not pipette by mouth. Avoid ingestion and contact with skin.

Serum

Blood should be collected from a patient and the serum separated as soon as possible. Soon after the blood is drawn, it should be allowed to clot, centrifuged, and the serum separated from the clot to a plastic tube (not glass) within 2 hours. It is recommended that specimen collection be carried out in accordance with NCCLS document M29-A2.

Avoid repeated freeze/thaw cycles.

Plasma

Whole blood is collected in sodium citrate, sodium EDTA and sodium fluoride anticoagulant. After collection, immediately centrifuge. In the U.S., follow NCCLS guideline H3-A2. Avoid repeated freeze/thaw cycles.

Serum or plasma may be stored refrigerated (2-8°C) for up to a week.¹ For long-term storage, keep at -20°C or below.

AUTOMATED ANALYZER APPLICATION

Suitable for two-reagent automated analyzers that use a multi-point calibration method.

PROCEDURE

Materials Supplied

Reagent 1 (R-1) Buffer Reagent 2 x 13.5 mL
Reagent 2 (R-2) Latex Suspension 2 x 5 mL

Materials Required But Not Supplied

Calibrators: **K-ASSAY®** Insulin Calibrator,
Cat. No. KAI-072C

Purified water

Clinical chemistry analyzer capable of accurately reading at 600 nm (main) and 800 nm (sub), accurately dispensing the required volumes, and maintaining 37°C.

Assay Procedure

Note: Allow all reagents and specimens to warm to room temperature. Mix all reagents gently before using.

An example of automated application:

Sample	12 μ L
↓	
• ←R1 (Buffer Reagent)	135 μ L
↓	37 °C, 5 min.
• ←R2 (Latex Suspension)	50 μ L
↓	37 °C, 5 min.

2 Point End at 600nm main, 800nm sub (if available)

Chemistry Parameters for Automatic Analyzer

INSTRUMENT	Roche / Hitachi 917
TEMPERATURE	37°C
TEST	(INS)
ASSAY CODE	(2 POINT END)(10) (19)(34)(0)(0)
WAVELENGTH	(800) (600)
SAMPLE VOLUME	(12.0) (0.0) (0)
REAGENT VOL (R1)	(135) (0)
REAGENT VOL (R2)	(0) (0)
REAGENT VOL (R3)	(50) (0)
REAGENT VOL (R4)	(0) (0)
ABS. LIMIT (SLOPE)	(32000) (INCREASE)
PROZONE LIMIT	(-32000) (34) (LOWER)
CALIB. TYPE	(SPLINE)
POINT	(6)
SPAN POINT	(6)
SD LIMIT	(999)
DUPLICATE LIMIT	(10000)
SENSITIVITY LIMIT	(0)
S1ABS RANGE	(-32000) (32000)
INSTRUMENT FACTOR	a=(1.0) b=(0.0)
UNIT	(μ U/mL)
STD.(1) Conc.-POS.	(0.0) - (1)
STD.(2) Conc.-POS.	(*2) - (2)
STD.(3) Conc.-POS.	(*3) - (3)
STD.(4) Conc.-POS.	(*4) - (4)
STD.(5) Conc.-POS.	(*5) - (5)
STD.(6) Conc.-POS.	(*6) - (6)

Use DI water for standard 1. Use Insulin calibrators 1-5 for standards 2-6.
*2-6: Input concentration of calibrators (using one decimal place [X.X]).

= User Defined

Parameters for other automated analyzers are available.

CALIBRATION

It is recommended that a multi-point calibration curve be made using the **K-ASSAY®** Insulin Calibrator (KAI-072C). It is recommended that the user determine calibration curve frequency as this depends on the instrument and type/number of other assays being performed. Initially, calibration should be performed each day.

QUALITY CONTROL

Normal and abnormal controls of known concentration should be included in each assay performed. These controls should fall within the ranges established by each lab for the particular lot of controls. The validity of the assay is in question if the value for the controls generated by the assay's calibration curve does not fall within the stated range. Recalibrate if the value determined for the controls falls outside the established recovery range.

RESULTS / CALCULATIONS

Insulin levels are determined using the prepared calibration curve.

LIMITATIONS OF PROCEDURE

The measurable range for insulin is 1.0 to 100 μ IU/mL. If the insulin value of a sample is greater than highest calibrator value, dilute 1 part sample with 3 parts isotonic saline and re-assay. Multiply results by 4 to compensate for the dilution.

PERFORMANCE

Specificity

When a sample with a known value is assayed, the result is within \pm 15% of the assigned value.

Precision

When a sample is assayed 5 times (within-run), the absorbance C.V. is \leq 10%.

(Within Run)

The following results were obtained on a Hitachi 917 analyzer with human serum:

	Sample I	Sample II	Sample III
N	10	10	10
Mean (μ IU/mL)	18.81	27.82	73.17
Std. Dev.	0.224	0.192	0.760
C.V. %	1.19	0.69	1.04

(Between Runs)

The following results were obtained on a Hitachi 917 analyzer with human serum:

	Sample IV	Sample V	Sample VI
N	10	10	10
Mean (μ IU/mL)	12.21	25.90	68.74
Std. Dev.	0.530	0.777	1.693
C.V. %	4.34	3.00	2.46

Accuracy / Correlation

A comparison of the **K-ASSAY**® Insulin and another company's Insulin EIA was performed with the following results:

Serum Sample

$$y = 0.8061x + 2.8674$$

$$r = 0.997$$

$$n = 32$$

x = another company's Insulin assay

$$y = \text{K-ASSAY}^{\circledR} \text{ Insulin Assay}$$

Plasma Sample

$$y = 0.8161x + 4.5552$$

$$r = 0.986$$

$$n = 47$$

x = another company's Insulin assay

$$y = \text{K-ASSAY}^{\circledR} \text{ Insulin Assay}$$

Assay Range

1.0 to 100 μ IU/mL (or value of highest calibration point)

Lower Limit of Detection

The analytical sensitivity is 1 μ IU/mL. This means that when saline and serum containing 1 μ IU/mL of insulin are tested 10 times, + 2.6 SD of the respective results do not overlap each other.

INTERFERENCE

No cross-reactivity with pro-insulin was observed. Hemoglobin, bile or rheumatoid factor did not interfere with the assay.

Bilirubin, Conjugated	No interference up to 19.9 mg/dL
Bilirubin, Unconjugated	No interference up to 19.3 mg/dL
Hemoglobin	No interference up to 450 mg/dL
Lipemia	No interference up to a formazin turbidity of 1,550
Rheumatoid Factor	No interference up to 4500 IU/L

PROZONE

No hook effect seen up to at least 1,000 μ IU/mL.











EXPECTED VALUES

The expected range for fasting insulin concentration has been reported to be up to 20 to 35 μ IU/mL (RIA).¹

REFERENCES

1. Jacobs DS, *et al.* Ed., Laboratory Test Handbook, 4th Edition. Lexicomp. p. 149 (1996).

LABELING SYMBOLS

	Catalog Number
	Expiration or "Use By" Date
	Lot Number
	Consult Package Insert for Instructions for Use
	For <i>In Vitro</i> Diagnostic Use
	CE Mark Registered
	For Prescription Use Only
	Temperature Limitation. Store between 2 and 8 degrees C
	Manufacturer
	Authorized Representative in the European Community

EU AUTHORIZED REPRESENTATIVE



Advena Ltd.

Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

ORDERING / PRICING / TECHNICAL INFORMATION



KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive
Seattle, WA 98168 USA
TEL: (206) 575-8068 / (800) 526-4925
FAX: (206) 575-8094