

K-ASSAY®

Insulin Calibrator

Lot AB123, Exp. 2025-02-07

Cat. No. KAI-072C

INTENDED USE

The **K-ASSAY®** Insulin Calibrator is a reagent for the calibration of the **K-ASSAY®** Insulin Assay. *FOR IN VITRO* DIAGNOSTIC USE.

SUMMARY

The calibrators in this kit contain known quantities of insulin. This kit is only to be used as a calibrator for the **K-ASSAY®** Insulin assay.

KIT COMPOSITION

Calibrator 1-5 5 x 0.5 mL, lyophilized

WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE. R only.

Do not mix or use calibrators from one test kit with those from a different lot number.

Do not use calibrators past their expiration date stated on each container label.

CALIBRATOR PREPARATION

1. Calibrators are lyophilized and require reconstitution. Add exactly 0.5 mL of purified water to each vial containing the lyophilized calibrator to dissolve it. Swirl gently to dissolve. The concentration of the dissolved calibrator is indicated on the vial label and may change from lot to lot.
2. Allow the vial to sit at room temperature (18-26°C) for 15 min. prior to use.
3. Mix before using by gently inverting the bottles.
4. Use the 5 reconstituted calibrators without dilution and purified water for the zero calibrator (0 µIU/mL) to make a calibration curve. Load the 6 calibrators from lowest to highest concentration on the analyzer.

STORAGE AND HANDLING

All calibrators should be stored refrigerated (2-8°C). Return all calibrators to 2-8°C promptly after use. Unopened calibrators can be used for 18 months from the date of manufacturing as indicated by the expiration date on the package and bottle labels. The opened calibrators are stable for two weeks when stored at 2-8°C.

CALIBRATOR STABILITY

Discard calibrators if they become contaminated. Evidence of cloudiness or particulate material in solution is cause to discard. Opened calibrators are stable for two weeks if stored at 2-8°C.

PROCEDURE**Materials Supplied**

Calibrator should be used as specified in the **K-ASSAY®** Insulin assay package insert.

Calibrator 1-5 5 x 0.5 mL, lyophilized

Materials Required But Not Supplied

K-ASSAY® Insulin assay, Cat. No. KAI-071

Pipette

Purified water

Details of Procedure

Note: Allow all reagents and specimens to come to room temperature (18-25°C). Mix all reagents gently before using.

K-ASSAY® Insulin Calibrators are assayed using the same procedure as the samples run in the test procedure. See package insert from the **K-ASSAY®** Insulin assay.

CALIBRATOR VALUES**Lot: AB123**

Calibrator 1 3.0 µIU/mL

Calibrator 2 8.2 µIU/mL











Calibrator 3 26.5 µIU/mL

Calibrator 4 51.1 µIU/mL

Calibrator 5 97.1 µIU/mL

The expected values for the **K-ASSAY®** Insulin Calibrator are continually being revised through ongoing quality assurance. As a result, the expected values may change from lot to lot. Please refer to the package insert for each lot for the exact calibrator values. Values are traceable to the Human Insulin for Immunoassay International Reference Preparation, NIBSC Code 66/304.

LABELING SYMBOLS

	Catalog Number
	Expiration or "Use By" Date
	Lot Number
	Consult Package Insert for Instructions for Use
	For <i>In Vitro</i> Diagnostic Use
	CE Mark Registered
	For Prescription Use Only
	Temperature Limitation. Store between 2 and 8 degrees C
	Manufacturer
	Authorized Representative in the European Community

EU AUTHORIZED REPRESENTATIVE



Advena Ltd.

Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

ORDERING / PRICING / TECHNICAL INFORMATION



KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive
Seattle, WA 98168 USA
TEL: (206) 575-8068 / (800) 526-4925
FAX: (206) 575-8094

Printed August 2023