

KAMIYA BIOMEDICAL COMPANY

ACE Inhibition Screening Kit

For the measurement of ACE inhibitory activity

Cat. No. KT-534

For Research Use Only. Not for Use in Diagnostic Procedures.

PRODUCT INFORMATION**ACE Inhibition Screening Kit**

Cat. No. KT-534

PRODUCT

The K-ASSAY® ACE Inhibition Screening Kit is for the measurement of ACE inhibitory activity.

BACKGROUND

Angiotensin-converting enzyme (ACE) is one of the key elements responsible for vasopressor action. ACE converts angiotensin I to angiotensin II, a potent vasopressor, in the renin-angiotensin system and contributes to increasing blood pressure by inactivating bradykinin, a strong antihypertensive peptide. Recently, various functional foods have received attention because of their inhibitory activity toward ACE.

ACE activity is conventionally determined by UV measurement of the hippuric acid produced from the synthetic substrate Hyppuryl-His-Leu. However, the assay process is complicated and requires organic solvent. In this kit, a safe and straightforward modified method has been developed.

The colorimetric detection system in the kit determines the amount of 3-hydroxybutyric acid (3HB) generated from 3-hydroxybutyryl-Gly-Gly-Gly by ACE. The kit is designed for 96-well microplate assays and is suitable for multiple sample measurements. No organic solvent extraction is required. The assay is safe, simple, and provides highly reproducible data.




Figure 1. Principle of ACE Inhibitory activity assay using the K-ASSAY® ACE Inhibition Screening Kit

COMPONENTS

- Substrate Buffer (1 mL), 2 vials
- Enzyme A, 2 vials
- Enzyme B, 2 vials
- Enzyme C, 2 vials
- Coenzyme, 2 vials
- Indicator solution (5 mL), 2 vials

MATERIALS OR EQUIPMENT REQUIRED BUT NOT PROVIDED

- Microplate reader (450 nm filter).
- 96-well microplate.
- 2 - 20 μ L, 20 - 200 μ L & 100 - 1,000 μ L pipettes
- Multi-channel pipette
- Incubator
- Disposable syringes (1 mL)

PREPARATION OF WORKING SOLUTION

A. Enzyme working solution

Dissolve Enzyme B in 2 mL of deionized water to prepare Enzyme B solution. Then add 1.5 mL of Enzyme B solution to Enzyme A to prepare Enzyme working solution.

Note: Enzyme A and B vials are capped under vacuum pressure. Add deionized water or solution through the rubber septum with a syringe, and then remove the septum.

Note: The Enzyme working solution is stable at -20°C for 2 weeks or in a refrigerator for 3 days.


B. Indicator working solution

Dissolve Enzyme C and Coenzyme in 3 mL of deionized water each. Add 2.8 mL of Enzyme C solution and 2.8 mL of Coenzyme solution to Indicator solution to prepare Indicator working solution.

Note: Enzyme C and Coenzyme vials are capped under vacuum pressure. Add deionized water through the rubber septum with a syringe, and then remove the septum.

Note: The Indicator working solution is stable at -20°C for 2 weeks or in a refrigerator for 3 days.

SELECTION OF PROTOCOL



PROTOCOL 1

Confirmation of the presence or absence of ACE inhibition

Sample Preparation

Prepare 100 μL of each sample.

Note: If the sample volume is $<100 \mu\text{L}$, please dilute the sample.

Note: If the sample is colored, prepare 150 μL of each sample.

Procedure for Colorless Sample

See Table 1 and Figure 2

- 1) Add 20 μL of sample solution to each sample well.
- 2) Add 20 μL of deionized water to each blank 1 well and 40 μL to each blank 2 well.
- 3) Add 20 μL of Substrate buffer to each well.
- 4) Add 20 μL of Enzyme working solution to each sample well and blank 1 well.

Note:

Enzyme working solution is easy to remain on the well wall. Please tap the plate to ensure that all solutions have been mixed completely.

Since the enzymatic reaction starts immediately after adding the Enzyme working solution, use a multi-channel pipette to minimize the well-to-well time lag.

- 5) Incubate at 37°C for 1 h.
- 6) Add 200 μL of Indicator working solution to each well.
- 7) Incubate at room temperature for 10 min.
- 8) Read the absorbance at 450 nm using a microplate reader.

Table 1 Amount of sample and reagent needed for each well.

	Sample	blank 1	blank 2
Sample solution	20 µl	-	-
Deionized water	-	20 µl	40 µl
Substrate buffer	20 µl	20 µl	20 µl
Enzyme working solution	20 µl	20 µl	-
Indicator working solution	200 µl	200 µl	200 µl

blank 1: positive control (no ACE inhibition)

blank 2: reagent blank

Figure 2 Example of arrangement on a 96-well microplate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1		Sample 8									
B	Sample 2		Sample 9									
C	Sample 3		Sample 10									
D	Sample 4		Sample 11									
E	Sample 5		Sample 12									
F	Sample 6		Sample 13									
G	Sample 7		Sample 14									
H	blank 1		blank 2									

Confirmation of the presence or absence of ACE inhibition

ACE inhibition can be calculated from the equation:

$$\text{ACE inhibition (\%)} = \frac{(\text{Ablank 1} - \text{Asample})}{(\text{Ablank 1} - \text{Ablank 2})} \times 100$$




Figure 3 Confirmation of the presence or absence of ACE inhibition.

Procedure for Colored Sample

See Table 2 and Figure 4

- 1) Add 20 µL of sample solution to each sample well and sample blank well.
- 2) Add 20 µL of deionized water to each blank 1 well, 40 µL to each blank 2 well and 240 µL to each sample blank well.
- 3) Add 20 µL of Substrate buffer to each sample well, blank 1 well and blank 2 well.
- 4) Add 20 µL of Enzyme working solution to each sample well and blank 1 well.

Note:

Enzyme working solution is easy to remain on the well wall. Please tap the plate to ensure that all solutions have been mixed completely.

Since the enzymatic reaction starts immediately after adding the Enzyme working solution, use a multi-channel pipette to minimize the well-to-well time lag.

- 5) Incubate at 37°C for 1 h.
- 6) Add 200 µL of Indicator working solution to each sample well, blank 1 well and blank 2 well.
- 7) Incubate at room temperature for 10 min.
- 8) Read the absorbance at 450 nm using a microplate reader.

Table 2 Amount of sample and reagent needed for each well.

	Sample	blank 1	blank 2	Sample blank
Sample solution	20 µl	-	-	20 µl
Deionized water	-	20 µl	40 µl	240 µl
Substrate buffer	20 µl	20 µl	20 µl	-
Enzyme working solution	20 µl	20 µl	-	-
Indicator working solution	200 µl	200 µl	200 µl	-

blank 1: positive control (no ACE inhibition)

blank 2: reagent blank

Figure 4 Example of arrangement on a 96-well microplate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	Sample 8	Sample blank 1	Sample blank 8								
B	Sample 2	Sample 9	Sample blank 2	Sample blank 9								
C	Sample 3	Sample 10	Sample blank 3	Sample blank 10								
D	Sample 4	Sample 11	Sample blank 4	Sample blank 11								
E	Sample 5	Sample 12	Sample blank 5	Sample blank 12								
F	Sample 6	Sample 13	Sample blank 6	Sample blank 13								
G	Sample 7	Sample 14	Sample blank 7	Sample blank 14								
H	blank 1	blank 2										

Confirmation of the presence or absence of ACE inhibition

ACE inhibition can be calculated from the equation:

$$\text{ACE inhibition (\%)} = [(A_{\text{blank 1}} - A_{\text{sample}} - A_{\text{sample blank}}) / (A_{\text{blank 1}} - A_{\text{blank 2}})] \times 100$$




Figure 5 Confirmation of the presence or absence of ACE inhibition

PROTOCOL 2

Determination of IC₅₀

Preparation of Sample solution

Dilute sample solution with deionized water.

Dilutions: 1 (no dilution), 1/5, 1/5², 1/5³, 1/5⁴, 1/5⁵, 1/5⁶




Figure 6. Preparation of sample solutions

Procedure for Colorless Sample

See Table 3 and Figure 7

- 1) Add 20 μL of sample solution to each sample well.
- 2) Add 20 μL of deionized water to each blank 1 well and 40 μL to each blank 2 well.
- 3) Add 20 μL of Substrate buffer to each well.
- 4) Add 20 μL of Enzyme working solution to each sample well and blank 1 well.

Note:

Enzyme working solution is easy to remain on the well wall. Please tap the plate to ensure that all solutions have been mixed completely.

Since the enzymatic reaction starts immediately after the addition of the Enzyme working solution, use a multi-channel pipette to minimize the well-to-well time lag.

- 5) Incubate at 37°C for 1 h.
- 6) Add 200 μL of Indicator working solution to each well.
- 7) Incubate at room temperature for 10 min.
- 8) Read the absorbance at 450 nm using a microplate reader.
- 9) ACE inhibition can be calculated from the equation:

$$\text{ACE inhibition (\%)} = \frac{[(\text{Ablank 1} - \text{Asample})/(\text{Ablank 1} - \text{Ablank 2})] \times 100}$$

Table 3 Amount of sample and reagent needed for each well.

	Sample	blank 1	blank 2
Sample solution	20 μl	-	-
Deionized water	-	20 μl	40 μl
Substrate buffer	20 μl	20 μl	20 μl
Enzyme working solution	20 μl	20 μl	-
Indicator working solution	200 μl	200 μl	200 μl

blank 1: positive control (no ACE inhibition)

blank 2: reagent blank


Figure 7 Example of arrangement on a 96-well microplate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	Sample 2										
B	Sample 1/5	Sample 2/5										
C	Sample 1/5 ²	Sample 2/5 ²										
D	Sample 1/5 ³	Sample 2/5 ³										
E	Sample 1/5 ⁴	Sample 2/5 ⁴										
F	Sample 1/5 ⁵	Sample 2/5 ⁵										
G	Sample 1/5 ⁶	Sample 2/5 ⁶										
H	blank 1	blank 2										

Determination of IC50 (50% inhibitory concentration)

- Prepare an inhibition curve using sample concentration for the x-axis and percentage inhibition of ACE for the y-axis. A typical inhibition curve is shown in Figure 8.
- Determine the concentration of the sample solution that gives 50% ACE inhibition as indicated in Figure 8.
- Because the total volume of the inhibition assay is 60 μL (first step of the assay), the original sample is diluted 3 times in the reaction. Therefore, the actual concentration of the sample at 50% inhibition is one-third of the concentration determined from the inhibition curve.

Figure 8 Inhibition curve



Procedure for Colored Sample

See Table 4 and Figure 9

- 1) Add 20 μ L of sample solution to each sample well and sample blank well.
- 2) Add 20 μ L of deionized water to each blank 1 well, 40 μ L to each blank 2 well and 240 μ L to each sample blank well.
- 3) Add 20 μ L of Substrate buffer to each sample well, blank 1 well and blank 2 well.
- 4) Add 20 μ L of Enzyme working solution to each sample well and blank 1 well.

Note:

Enzyme working solution is easy to remain on the well wall. Please tap the plate to ensure that all solutions have been mixed completely.

Since the enzymatic reaction starts immediately after adding the Enzyme working solution, use a multi-channel pipette to minimize the well-to-well time lag.

- 5) Incubate at 37°C for 1 h.
- 6) Add 200 μ L of Indicator working solution to each sample well, blank 1 well and blank 2 well.
- 7) Incubate at room temperature for 10 min.
- 8) Read the absorbance at 450 nm using a microplate reader.
- 9) ACE inhibition can be calculated from the equation.

$$\text{ACE inhibition (\%)} = \frac{(\text{Ablank 1} - \text{Asample} - \text{Asample blank})}{(\text{Ablank 1} - \text{Ablank 2})} \times 100$$

Table 4 Amount of sample and reagent needed for each well.

	Sample	blank 1	blank 2	Sample blank
Sample solution	20 μ l	-	-	20 μ l
Deionized water	-	20 μ l	40 μ l	240 μ l
Substrate buffer	20 μ l	20 μ l	20 μ l	-
Enzyme working solution	20 μ l	20 μ l	-	-
Indicator working solution	200 μ l	200 μ l	200 μ l	-

blank 1: positive control (no ACE inhibition)

blank 2: reagent blank

Figure 9 Example of arrangement on a 96-well microplate


	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	Sample 2	Sample blank 1	Sample blank 2								
B	Sample 1/5	Sample 2/5	Sample blank 1/5	Sample blank 2/5								
C	Sample 1/5 ²	Sample 2/5 ²	Sample blank 1/5 ²	Sample blank 2/5 ²								
D	Sample 1/5 ³	Sample 2/5 ³	Sample blank 1/5 ³	Sample blank 2/5 ³								
E	Sample 1/5 ⁴	Sample 2/5 ⁴	Sample blank 1/5 ⁴	Sample blank 2/5 ⁴								
F	Sample 1/5 ⁵	Sample 2/5 ⁵	Sample blank 1/5 ⁵	Sample blank 2/5 ⁵								
G	Sample 1/5 ⁶	Sample 2/5 ⁶	Sample blank 1/5 ⁶	Sample blank 2/5 ⁶								
H	blank 1	blank 2										

Determination of IC₅₀ (50% inhibitory concentration)

- Prepare an inhibition curve using sample concentration for the x-axis and percentage inhibition of ACE for the y-axis. A typical inhibition curve is shown in Figure 10.

- Determine the concentration of the sample solution that gives 50% ACE inhibition as indicated in Figure 10.
- Because the total volume of the inhibition assay is 60 μL (first step of the assay), the original sample is diluted 3 times in the reaction. Therefore, the actual concentration of the sample at 50% inhibition is one-third of the concentration determined from the inhibition curve.

Figure 10 Inhibition curve



STORAGE

Store the kit at 4°C.

PRECAUTION

- Several kit components are in glass vials. Please handle with care.
- Multiple measurements (triplicates of each sample) are recommended to obtain accurate data.
- If the water solubility of the sample is low, use dimethylsulfoxide or ethanol to dissolve. Then, dilute the solution with an appropriate buffer. The final concentration of organic solvent should be <1%.
- If the sample solution is acidic, adjust the pH to ≥ 5 before use for measurement.
- Ascorbic acid may interfere with the assay. The concentration of ascorbic acid in the sample solution should be <0.01% w/v. If the sample solution contains insoluble materials, remove it by centrifugation or filtration before use for measurement.

**FOR RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES**

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168
Tel: (206) 575-8068 Fax: (206) 575-8094
Email: LifeScience@k-assay.com
www.k-assay.com